

Cell culture and adipocyte differentiation

 Marcelo A. Mori  Andréa Livia Rocha

Updated date: Dec 21, 2020



An abbreviated version of this protocol was published in Science Advances in Dec 2020

Enoxacin induces oxidative metabolism and mitigates obesity by regulating adipose tissue miRNA expression

DOI: 10.1126/sciadv.abc6250

Detailed protocol

Materials and Reagents

1. Murine 3T3-F442A, C3H10T1/2, white subcutaneous (9W), or brown (9B) preadipocytes
2. 20 nM insulin (Sigma-Aldrich, catalog number: I6634)
3. 1 nM triiodothyronine (T3) (Sigma-Aldrich, catalog number: T2877)
4. 0.5 mM 3-Isobutyl-1-methylxanthine (Sigma-Aldrich, catalog number: I5879)
5. 1 μ M dexamethasone (Sigma-Aldrich, catalog number: D4902)
6. 0.125 mM indomethacin (Sigma-Aldrich, catalog number: I7378)
7. 2.8 μ M rosiglitazone (Sigma-Aldrich, catalog number: R2408)
8. Dulbecco's modified Eagle's medium (high glucose) (Life Technologies, Invitrogen™, catalog number: 11965)
9. Fetal bovine serum (FBS) (Thermo Fisher Scientific, catalog number: 12657-029)
10. Pen/Strep (Thermo Fisher Scientific, catalog number: 15140122)

Equipment

1. 37 °C, 5% CO₂ incubator
2. Water bath

Procedure

Notes:

Pre-warm DMEM media (see Recipes) to 37 °C prior to the experiment

A. Seed the cells

1. Seed 1×10^5 preadipocytes in a 12 well plate and cultivate in DMEM media (see Recipes) until 90% confluence
2. Change the media every 2 days

B. Adipocyte differentiation (pro-browning stimulus)

- Day 0: 1. Once the cells reach 90% confluence, change the media and cultivate for 2 days
- Day 2 2. Replace the media for Differentiation media day 2 (see Recipes) and cultivate for 2 days
- Day 4 3. Remove the media and replace with Differentiation media day 4 (see Recipes)
- Day 6 4. After 2 days, replace the media with Differentiation media day 6 (see Recipes)
- Day 8 5. Cells are fully differentiated and ready to use

Recipes

1. DMEM media (supplemented with 10% FBS and 1% of Pen/Strep)
2. Differentiation media day 2: DMEM media further supplemented with Insulin, Dexamethasone, IBMX, T3, Indomethacin, Rosiglitazone (final concentrations in Material and Reagents)
3. Differentiation media day 4: DMEM media further supplemented with Insulin, T3 and Rosiglitazone (final concentrations in Material and Reagents)
4. Differentiation media day 6: DMEM media further supplemented with Insulin and Rosiglitazone (final concentrations in Material and Reagents)

Acknowledgments

We are grateful for the technical assistance of E. E. Saviani and for the supported by grants from FAPESP (2017/01184-9, 2017/07975-8, 2016/12294-7, 2010/52557-0, 2016/24163-4, 2016/23142-3, 2018/15313-8, and 2015/01316-7), CNPq (474397/2011-4 and 305069/2015-2), CAPES (1431744/2014-2016), and the Ecole Polytechnique Fédérale de Lausanne. INFABIC is co-funded by FAPESP and CNPq.

Representative images




Figure 1. Representative microscopic images of 9W preadipocytes before and after differentiation. Image of 9W preadipocytes at day 0 - 100% confluence and at day 8 - fully differentiated adipocytes.

How to cite: (Readers should cite both the Bio-protocol preprint and the original research article where this protocol was used)

1. Mori, M. A. and Rocha, A. L.(2020). Cell culture and adipocyte differentiation. Bio-protocol Preprint. [bio-protocol.org/prep712](https://doi.org/10.21956/bio-protocol.d712).
2. Rocha, A. L., Lima, T. I. D., Souza, G. P. D., Corrêa, R. O., Ferrucci, D. L., Rodrigues, B., Lopes-Ramos, C., Nilsson, D., Knittel, T. L., Castro, P. R., Fernandes, M. F., Martins, F. D. S., Parmigiani, R. B., Silveira, L. R., Carvalho, H. F., Auwerx, J., Vinolo, M. A. R., Boucher, J. and Mori, M. A.(2020). Enoxacin induces oxidative metabolism and mitigates obesity by regulating adipose tissue miRNA expression. Science Advances 6(49). DOI:

Copyright: Content may be subjected to copyright.